

# Initial F MTP™: Single use sterile kits for 1<sup>st</sup> metatarso-phalangeal arthrodesis

## Initial F MTP™ - Kits references



### Indications:

The implant of the Initial F MTP™ is intended for arthrodeses, fractures and osteotomies fixation and revision surgeries of the foot in adults.

### Contra-indications:

- Serious vascular deterioration, bone devitalization,
- Pregnancy,
- Acute or chronic local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency affecting the concerned area.
- Insufficient bone quality preventing a good fixation of the implants into the bone,
- Muscular deficit, neurological deficiency or behavioral disorders, which could submit the implant to abnormal mechanical strains.
- Allergy to one of the materials used or sensitivity to foreign bodies.
- Serious problems of non-compliance, mental or neurological disorders, failure to follow post-operative care recommendations.
- Unstable physical and/or mental condition.

**SLT2.8Lxx**  
Ø2.8 mm locking screw  
Green anodized


Pin Ø1.6 L100 mm

Pin Ø1.2 L70 mm  
x2

T8 prehensor  
screwdriver

Length gauge

Ø2.0 mm quick  
coupling drill bit  
L 125 mm



**RLT2.8Lxx**  
Ø2.8 mm non locking screw  
Yellow anodized

Handle  
for guide gauge

Ø2.0 mm threaded  
guide gauge

# Initial F<sup>TM</sup> MTP - Kits references

## > Kits content - Implants



1<sup>st</sup> MTP arthrodesis plate  
Size 1



1<sup>st</sup> MTP arthrodesis plate  
Size 2

KIT-MD1D or KIT-MD1G* 1st MTP arthrodesis kit - Left or Right - Size 1		
Ref.	Description	KIT-MD1D or KIT-MD1G
FMTxD1	1 <sup>st</sup> MTP arthrodesis plate - Size 1 - Right or Left	1
RLT2.8L16	Non locking screw - Ø2.8 mm - L 16 mm	1
RLT2.8L18	Non locking screw - Ø2.8 mm - L 18 mm	1
RLT2.8L22	Non locking screw - Ø2.8 mm - L 22 mm	1
RLT2.8L26	Non locking screw - Ø2.8 mm - L 26 mm	1
SLT2.8L14	Locking screw - Ø2.8 mm - L14 mm	1
SLT2.8L16	Locking screw - Ø2.8 mm - L16 mm	2
SLT2.8L18	Locking screw - Ø2.8 mm - L18 mm	2
SLT2.8L20	Locking screw - Ø2.8 mm - L20 mm	1

KIT-MD2D or KIT-MD2G* 1st MTP arthrodesis kit - Left or Right - Size 2		
Ref.	Description	KIT-MD2D or KIT-MD2G
FMTxD2	1 <sup>st</sup> MTP arthrodesis plate - Size 2 - Right or Left	1
RLT2.8L16	Non locking screw - Ø2.8 mm - L 16 mm	1
RLT2.8L18	Non locking screw - Ø2.8 mm - L 18 mm	1
RLT2.8L22	Non locking screw - Ø2.8 mm - L 22 mm	1
RLT2.8L26	Non locking screw - Ø2.8 mm - L 26 mm	1
SLT2.8L14	Locking screw - Ø2.8 mm - L14 mm	1
SLT2.8L16	Locking screw - Ø2.8 mm - L16 mm	2
SLT2.8L18	Locking screw - Ø2.8 mm - L18 mm	2
SLT2.8L20	Locking screw - Ø2.8 mm - L20 mm	1

## > Additional instrumentation kits


### SINGLE USE CONVEX AND CONCAVE REAMERS - STERILE PACKAGING

Ref.	Description	Content
KIT-MI16	Ø16 mm reamers kit for 1st MTP arthrodesis	Ø16 mm Initial convex reamer Ø16 mm Initial concave reamer Pin Ø1.6 L100 mm x 2
KIT-MI19	Ø19 mm reamers kit for 1st MTP arthrodesis	Ø19 mm Initial convex reamer Ø19 mm Initial concave reamer Pin Ø1.6 L100 mm x 2
KIT-MI22	Ø22 mm reamers kit for 1st MTP arthrodesis	Ø22 mm Initial convex reamer Ø22 mm Initial concave reamer Pin Ø1.6 L100 mm x 2

### INITIAL F<sup>TM</sup> - MTP TEMPLATES

Ref.	Description
ANC808*	1 <sup>st</sup> MTP arthrodesis plates template
ANC846	1 <sup>st</sup> MTP arthrodesis reamers template

#### > Reamers



#### > Template plate - ANC808



\* Seules les tailles 1 et 2 sont disponibles en version Initial F MTP

\* Each kit includes one plate (left or right).  
NB: Additional sterile screws can be ordered separated (Supplemental sterile screw caddy).



NEWCLIP TECHNICS  
PA de la Lande Saint Martin - 45 rue des Garotières  
44115 Haute Goulaine (France)  
Tél. : +33 (0)2 28 21 37 12 - Fax : +33 (0)2 40 63 68 37  
Mail : orders@newcliptechnics.com  
www.newcliptechnics.com

NEWCLIP USA  
642 Larkfield Center  
Santa Rosa CA 95403, USA  
Phone : + 1 707 230 5078  
Mail : customerservice@newclipusa.com  
www.newclipusa.com

The information presented in this brochure is intended to demonstrate a Newclip Technics product. Always refer to the package insert, product label and/or user instructions before using any Newclip Technics product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Newclip Technics representative if you have questions about the availability of Newclip Technics products in your area.

